

## Purification of Plasmids for Sequencing

*The following protocol for preparing good quality DNA template for sequencing is recommended by Applied Biosystems and can be found online in their Automatic DNA Sequencing Guide at:*

<http://www.pebi-docs.com/pebi-docs/04305080.pdf>

When purifying recombinant plasmids in bacteria, plate out the transformants to obtain isolated colonies. Select a single colony and re-streak it on a plate. Select an isolated colony from that plate to obtain plasmids with the desired insert.

The optimal procedure for preparing a particular plasmid depends on the particular bacterial strain and the yield of each construct. Good sequencing data has been obtained from the following methods:

- \* ABI PRISM™ Plasmid Miniprep Kit (P/N 402790 or 402791)
- \* Cesium chloride (CsCl) banding
- \* Modified alkaline lysis/PEG precipitation method (see below)
- \* PureGene DNA Isolation Kit (Gentra Systems, Inc., P/N D-5500A)
- \* QIAGEN Plasmid Kits (<http://www.qiagen.com>):
  - Mini (P/N 12123, 25 reactions; 12125, 100 reactions)
  - Midi (P/N 12143, 25 reactions; 12144, 50 reactions; 12145, 100 reactions)
  - Maxi (P/N 12162, 10 reactions; 12163, 25 reactions; 12165, 100 reactions)

### Modified Alkaline Lysis/PEG Method

Reagents required:

- \* Chloroform (or chloroform:isoamyl alcohol (24:1)).
- \* Deionized water
- \* Ethanol, 70%
- \* GET buffer (50 mM glucose, 10 mM EDTA, 25 mM Tris, pH 8.0)
- \* Isopropanol, 100% (anhydrous)
- \* PEG 8000, 13% (sterilized by autoclaving, rather than by filtration)
- \* Potassium acetate, 3 M, pH 4.8
- \* RNase A (DNase-free), 10 mg/mL
- \* Sodium chloride (NaCl), 4 M
- \* Sodium hydroxide (NaOH), 0.2 N, with 1% SDS (freshly made)

**Note** To minimize shearing of contaminating chromosomal DNA, **do not use a vortex during this procedure.**

1. Pellet 1.5-mL aliquots of culture for 1 minute in a microcentrifuge at maximum speed.

*Note:* A total culture volume of 4.5 mL can be spun down per tube without changing volumes in the procedure. This allows you to achieve a threefold increase in yield while eliminating the need for extra tubes and additional handling.

2. Remove the supernatant by aspiration.

3. Resuspend the bacterial pellet in 200  $\mu$ L of GET buffer by pipetting up and down.

4. Add 300  $\mu$ L of freshly prepared 0.2 N NaOH/1% SDS. Mix the contents of the tube by inversion. Incubate on ice for 5 minutes.

5. Neutralize the solution by adding 300  $\mu$ L of 3.0 M potassium acetate, pH 4.8. Mix by inverting the tube. Incubate on ice for 5 minutes.

6. Remove cellular debris by spinning in a microcentrifuge at maximum speed for 10 minutes at room temperature. Transfer the supernatant to a clean tube.

7. Add RNase A (DNase-free) to a final concentration of 20  $\mu$ g/mL. Incubate the tube at 37 °C for 20 minutes.

*\*See Pat's Note, below.*

8. Extract the supernatant twice with chloroform:

a. Add 400  $\mu$ L of chloroform.

b. Mix the layers by inversion for 30 seconds.

c. Centrifuge the tube for 1 minute to separate the phases.

d. Transfer the upper aqueous phase to a clean tube.

9. Add an equal volume of 100% isopropanol. Mix the contents of the tube by inversion.

10. Spin the tube in a microcentrifuge at maximum speed for 10 minutes at room temperature.

11. Remove the isopropanol completely by aspiration.

12. Wash the DNA pellet with 500  $\mu$ L of 70% ethanol. Dry under vacuum for 3 minutes.

13. Dissolve the pellet in 32  $\mu$ L of deionized water.

14. Add 8.0  $\mu$ L of 4 M NaCl, then 40  $\mu$ L of autoclaved 13% PEG 8000.

15. Mix thoroughly, then leave the sample on ice for 20 minutes.

16. Pellet the plasmid DNA by spinning in a microcentrifuge at maximum speed for 15 minutes at 2–6 °C.

17. Carefully remove the supernatant. Rinse the pellet with 500  $\mu$ L of 70% ethanol.

18. Resuspend the pellet in 20  $\mu$ L of deionized water. Store at –15 to –25 °C.

*Pat's Note: After step 7 the procedure may be modified to include a phenol extraction:*

*7a. Add an equal volume of phenol:chloroform:isoamyl alcohol (25:24:1) to the tubes, and mix by inversion until an emulsion forms.*

*7b. Centrifuge the mixture for 1 min. in a microfuge at room temperature. Use a pipette to transfer the aqueous phase to a fresh tube. Discard the interface and organic phase,*

*7c. Repeat this extraction. (Continue on to step 8.)*